**OPED GmbH**  
Medizinpark 1  
83626 Valley/Oberlaindern  
Germany  
Fon +49 (0) 80 24/60 81 8-210  
Fax +49 (0) 80 24/60 81 8-299  
mail@oped.de  
www.oped.de

**OPED AG**  
Hinterbergstrasse 26  
6330 Cham  
Switzerland  
Fon +41 (0) 41/748 53 35  
Fax +41 (0) 41/748 53 40  
contact@oped.ch  
www.oped.ch

**OPED UK Ltd**  
Unit 2 Glenmore Business Centre  
Waller Road  
Hopton Industrial Estate  
Devizes  
SN10 2EQ  
United Kingdom  
Phone +44 (0) 1380 722177  
Fax +44 (0) 1380 710095  
enquiries@oped.biz  
www.opeduk.co.uk

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Ihr OPED-Team

## ADVAGOshoulder

Intended use: orthosis to immobilise the shoulder joint in adjustable positions on two levels.

### Safety Advice

- ⚠ The product is single use and should only be used by one patient. 
- Should it be used for the treatment of more than one patient, the manufacturer's product liability shall lapse in accordance with the Medical Devices Act.
- ⚠ Never apply the product to open wounds and only use under prior medical instruction. 
- ⚠ The patient may only adjust it in consultation with his/her treating doctor.
- ⚠ Should you experience excessive pain, swelling, numbness or sensitivity changes, unusual reactions or an uncomfortable sensation while wear-

ing the ADVAGOshoulder, please contact your doctor or orthopaedic technician immediately.

- ⚠ This orthosis was developed for the indications stated and may be damaged through unusual and non-intended use or through activities such as sport, falls or other accidents, as well as liquids.
- ⚠ Please consult your treating doctor before using it in conjunction with another medical device.
- ⚠ Please observe the road transport law of your country.
- ⚠ Note: Not made with natural rubber latex. 

### How to Use

- Move sideways towards the ADVAGOshoulder and place the arm in the arm support (b). Fasten velcro strap (b). (Fig. 1-2)
- Place the hip belt (c) around the stomach and fasten. Tighten hip belt. The base (d) of the arm support should rest on the pelvic bone. (Fig. 3-4)
- Slip the healthy arm into the shoulder strap (e) and place it over the healthy shoulder. (Fig. 5-6)
- Insert hook (f) into the arm support. Pull the open belt clip of the ABC belt clip system (g) upwards until the breast strip is taut, then fasten the clip. (h) must always rest on the breastbone. (Fig. 7-8, 10)

The length of the shoulder strap can be adjusted using the green belt clasp in the shoulder pad. (Fig. 9)

- Arm support length adjustment: Press the button (h) down and adjust to the desired setting.
- Hand support length adjustment: Press the clips (i) on the bottom together and slide them into the desired position. The hand support is detachable.

### Cleaning and Care

**Washing:** Wash strap and pad in the washing machine at 40°C.  
Clean plastic parts with a damp cloth and with plastic-friendly cleaning agents only.

**Drying:** leave to dry naturally

### Product Characteristics

**Sizes:** Medium, large, X large, for left and right respectively

### Indications

- Rupture of the rotator cuff
- Tendon ruptures (supraspinatus, biceps)
- Endoprostheses
- AC joint separation
- Fractures of the greater tubercle
- Dislocations

### Contraindications

None known

### Technical data/parameters

- **Material composition:**  
Arm support: PP, small parts: PA, POM strap, pad: PP, PS, cotton
- **Dimensions:**  
(without strap): 44cm x 22cm x 21cm
- **Weight:** 1,050 g

### Transport conditions/ storage conditions/disposal:

**Storage:** Store in a dry place and out of direct sunlight

### Guarantee

The guarantee is in line with statutory prescriptions.

### Patient information

#### Possible adjustments:

Only in consultation with your doctor/therapist.

#### Internal/external rotation:

Pull the lever (j) inwards and slide it into the desired position to set the angle adjustment. Repeat this procedure for the other direction.

#### Adduction/abduction:

Open the flap (k) and adjust to the desired angle. To set the adjustment, press the flap down again. Please note that the chest strap also needs to be adjusted.

#### Showering/hygiene:

You can shower/be in water while wearing the product, subject to the approval of your treating doctor.

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